


GOAT ANTI-IL6R / CD126 (AA249-263) ANTIBODY

SKU: EB11139



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	IL6R OTTHUMP00000217870 MGC104991 membrane glycoprotein 80 interleukin-6 receptor subunit alpha interleukin 6 receptor alpha subunit interleukin 6 receptor IL-6R-alpha IL-6RA IL6RA IL-6R-1 IL-6R subunit alpha IL-6R 1 IL-6 receptor subunit alpha gp80 CD126 antigen CD126
Accession ID	NP_000556.1; NP_852004.1
Blocking Peptide	EBP11139
Immunogen	Peptide with sequence C-YRLRFELRYAERSK, from the internal region of the protein sequence according to NP_000556.1; NP_852004.1.
Product Comments	This antibody is expected to recognize all reported isoforms (NP_000556.1; NP_852004.1). The sequence of the immunizing peptide corresponds with a part of the extracellular domain.
Peptide Sequence	C-YRLRFELRYAERSK
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Pig
Reactive Species	Human, Pig
Human Gene ID	3570
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 45kDa band observed in Human and Pig Liver lysates (calculated MW of 40.2kDa according to NP_852004.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

