


Telephone: (650) 697-3600

GOAT ANTI-INHIBIN ALPHA CHAIN ANTIBODY

SKU: EB11943



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing. Storage

Synonym /

Alias inhibin alpha chain|OTTHUMP00000065425|A-inhibin subunit|inhibin, alpha|INHA

Names

Immunofluorescence: Strong expression of the protein seen in HeLa and U2OS cells.

Usage Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of HeLa **Summary**

cells. Recommended concentration: 10ug/ml.


Accession

NP 002182.1

Blocking EBP11943 Peptide









Telephone: (650) 697-3600

Peptide with sequence C-RPLHVRTTSDGGYS, from the C Terminus of the protein sequence according to **Immunogen**

NP_002182.1.

Peptide

C-RPLHVRTTSDGGYS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Refrigerated Instructions

Predicted Species

Human, Cow

Reactive

Human **Species**

Human

3623 **Gene ID**

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png

ELISA

Detection

Antibody detection limit dilution 1:128000.

Limit


Western Approx. 40kDa band observed in lysates of cell lines HepG2 and HeLa (calculated MW of 39.7kDa according to Blot NP_002182.1). Recommended concentration: 0.5-1µg/ml. Primary incubation 1 hour at room temperature.

Application


Type

Pep-ELISA, WB, IF, FC

GALLERY IMAGES











Telephone: (650) 697-3600



