



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-INO80 ANTIBODY**

**SKU:** EB09691

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

homolog of yeast INO80|INO80 complex subunit A|INO80 complex homolog 1|hINO80|INOC1|INO80A|INO80 **Alias** 

homolog (S. cerevisiae)|INO80 **Names** 

Accession

NP 060023.1

**Blocking** 

EBP09691 **Peptide** 

Peptide with sequence C-TQEDYDSNHFKAQ, from the internal region of the protein sequence according to Immunogen

NP 060023.1.

**Peptide** C-TQEDYDSNHFKAQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Rat Dog, Cow **Species** 

Human

54617 **Gene ID** 

**Rat Gene ID** 296084

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png Grade

**ELISA** 

**Detection** Antibody detection limit dilution 1:4000.

Limit

**Application** Pep-ELISA **Type** 

> Goat Anti-INO80 Antibody https://everestbiotech.com/products/goat-anti-ino80-antibody/