

GOAT ANTI-IQGAP1 ANTIBODY

SKU: EB09580



250kDa
150kDa
100kDa
75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	RasGAP-like with IQ motifs p195 SAR1 KIAA0051 HUMORFA01 IQ motif containing GTPase activating protein 1 IQGAP1
Accession ID	NP_003861.1
Blocking Peptide	EBP09580
Immunogen	Peptide with sequence C-KRQSGQTDPLQKEE, from the internal region of the protein sequence according to NP_003861.1.
Peptide Sequence	C-KRQSGQTDPLQKEE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	8826
Mouse Gene ID	29875
Rat Gene ID	361598
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 200kDa band observed in Human Kidney lysates (calculated MW of 189kDa according to NP_003861.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa