

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-IQGAP1 ANTIBODY**

**SKU:** EB09580







250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

RasGAP-like with IQ motifs|p195|SAR1|KIAA0051|HUMORFA01|IQ motif containing GTPase activating protein

1|IQGAP1 Names

**Accession** 

NP\_003861.1 ID

**Blocking** 

EBP09580 **Peptide** 

Peptide with sequence C-KRQSGQTDPLQKEE, from the internal region of the protein sequence according to **Immunogen** 

NP 003861.1.

**Peptide** Sequence

C-KRQSGQTDPLQKEE

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat **Species** 

Reactive

Human **Species** 

Human

8826 Gene ID

Mouse

29875 **Gene ID** 

**Rat Gene ID** 361598

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**Grade ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

Western Approx 200kDa band observed in Human Kidney lysates (calculated MW of 189kDa according to NP 003861.1).

Blot Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB

**Type** 



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GALLERY IMAGES**

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa