




GOAT ANTI-IQGAP2 ANTIBODY

SKU: EB09511



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias

IQ motif containing GTPase activating protein 2|IQGAP2

Names

Accession

NP_006624.2

Blocking

Peptide

EBP09511

Immunogen

Peptide with sequence C-NTRRSIKLDGKGEP, from the internal region (near C Terminus) of the protein sequence according to NP_006624.2.

Peptide

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human, Mouse, Dog, Cow

Species Reactive

Species Human

Human

Gene ID

10788

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results

 $In \ paraffin \ embedded \ Human \ Test is \ shows \ membranous \ and \ cytoplasm \ staining \ of \ primary \ spermatocytes.$

Recommended concentration, 2-6µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Application

Pep-ELISA, IHC

Type

GALLERY IMAGES

