

GOAT ANTI-IRE1 / ERN1 ANTIBODY

SKU: EB05645



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	protein kinase/endoribonuclease inositol-requiring 1 ER to nucleus signalling 1 MGC163279 MGC163277 IRE1P IRE1 FLJ30999 endoplasmic reticulum to nucleus signalling 1 ERN1
Accession ID	NP_001424.3
Blocking Peptide	EBP05645
Immunogen	Peptide with sequence C-PPEPQPPVTPDAL, from the C Terminus of the protein sequence according to NP_001424.3.
Peptide Sequence	C-PPEPQPPVTPDAL
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	2081
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Preliminary experiments in Human Muscle, Pancreas and Placenta lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA