


GOAT ANTI-IRF4 ANTIBODY

SKU: EB11330



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	IRF4 MUM1 multiple myeloma oncogene 1 LSIRF interferon regulatory factor 4
Accession ID	NP_002451.2
Blocking Peptide	EBP11330
Immunogen	Peptide with sequence C-HISNPEDYHRSIR, from the C Terminus of the protein sequence according to NP_002451.2.
Peptide Sequence	C-HISNPEDYHRSIR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human, Mouse
Human Gene ID	3662
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 48kDa band observed in lysates of cell lines Daudi and Jurkat and approx,48-50kDa in lysates of cell line NIH3T3 (calculated MW of 51.8kDa according to Human NP_002451.2 and 51.4kDa according to Mouse NP_001334437.1). Recommended concentration: 0.5-2µg/ml. Primary incubation 1 hour at room temperature.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

