

GOAT ANTI-JAM2 / JAMB / CD322 ANTIBODY

SKU: EB09490



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	JAM2 junctional adhesion molecule 2 C21orf43 CD322 JAM-B JAMB PRO245 VE-JAM VEJAM JAM-IT VE-JAM OTTHUMP0000096100 junctional adhesion molecule B vascular endothelial junction-associated molecule
Accession ID	NP_067042.1
Blocking Peptide	EBP09490
Immunogen	Peptide with sequence C-QRKGYSKETSFQ, from the internal region of the protein sequence according to NP_067042.1.
Peptide Sequence	C-QRKGYSKETSFQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	58494
Mouse Gene ID	67374
Rat Gene ID	619374
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Approx 35kDa band observed in Human Heart lysates (calculated MW of 33.2kDa according to NP_067042.1). Recommended concentration: 0.5-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

