

# GOAT ANTI-KAISO ANTIBODY

**SKU:** EB05263



## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	WUGSC:H_DJ525N14.1 KAISO Kaiso ZNF-kaiso
<b>Accession ID</b>	NP_006768
<b>Blocking Peptide</b>	EBP05263
<b>Immunogen</b>	Peptide with sequence C-DGSTEFEFIIPESY, from the C Terminus of the protein sequence according to NP_006768.
<b>Peptide Sequence</b>	C-DGSTEFEFIIPESY
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human, Mouse
<b>Reactive Species</b>	Mouse
<b>Human Gene ID</b>	10009
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>Western Blot</b>	Western blot: Approx 70kDa band seen in Mouse Spleen extracts. Recommended for use at 1-3µg/ml.