


GOAT ANTI-KATNAL2 ANTIBODY

SKU: EB12051



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	p60 katanin-like 2 katanin p60 ATPase-containing subunit A-like 2 katanin p60 subunit A-like 2 KATNAL2
Accession ID	NP_112593.2
Blocking Peptide	EBP12051
Immunogen	Peptide with sequence C-TGDTKSLNKEHPNQE, from the internal region of the protein sequence according to NP_112593.2.
Peptide Sequence	C-TGDTKSLNKEHPNQE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	83473
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 55kDa band observed in Human Testis lysates (calculated MW of 52.8kDa according to NP_112593.2). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

