

GOAT ANTI-KCC4 / SLC12A7 ANTIBODY

SKU: EB08089



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	potassium/chloride transporter KCC4 KCC4 DKFZP434F076 solute carrier family 12 (potassium/chloride transporters), member 7 SLC12A7
Accession ID	NP_006589.2
Blocking Peptide	EBP08089
Immunogen	Peptide with sequence C-QLIHDRNTASHT, from the internal region of the protein sequence according to NP_006589.2.
Peptide Sequence	C-QLIHDRNTASHT
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	10723
Mouse Gene ID	20499
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Preliminary experiments in Human Heart and Kidney lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

