

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-KCC4 / SLC12A7 ANTIBODY

SKU: EB08089

×

SPECIFICATIONS

Unit Size

100 μg

Synonym / Alias

 $potassium/chloride\ transporter\ KCC4|KCC4|DKFZP434F076|solute\ carrier\ family\ 12\ (potassium/chloride\ potassium/chloride\ potassium/chloride$

Names transporters), member 7|SLC12A7

Accession

ID

NP_006589.2

Blocking Peptide

EBP08089

Immunogen

Peptide with sequence C-QLIHDRNTASHT, from the internal region of the protein sequence according to

NP 006589.2.

Peptide

Sequence

C-QLIHDRNTASHT

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Rat

Human

Human Gene ID 10723

Ochic ib

Mouse 20499

Gene ID

iene ib

Product

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws. com/assets/products/image/aspiring_medium.png$

ELISA

Detection Limit Antibody detection limit dilution 1:2000.

Western

Blot

Preliminary experiments in Human Heart and Kidney lysates gave no specific signal but low background (at antibody concentration up to $1\mu g/ml$). We would appreciate any feedback from people in the field - have any

results been reported with other antibodies/lysates?

Application

Type

Pep-ELISA