


GOAT ANTI-KEAP1 (AA381-394) ANTIBODY

SKU: EB10962



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	kelch-like protein 19 cytosolic inhibitor of Nrf2 MGC4407 MGC20887 MGC1114 MGC10630 KLHL19 KIAA0132 INrf2 kelch-like ECH-associated protein 1 KEAP1
Usage Summary	Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: 10ug/ml.
Accession ID	NP_036421.2
Blocking Peptide	EBP10962
Immunogen	Peptide with sequence C-NNSPDGNTDSSALD, from the internal region of the protein sequence according to NP_036421.2.
Product Comments	Reported variants represent identical protein: NP_036421.2, NP_987096.1
Peptide Sequence	C-NNSPDGNTDSSALD

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	9817
Mouse Gene ID	50868
Rat Gene ID	117519
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA, FC

GALLERY IMAGES

