

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-KI-67 / MKI67 ANTIBODY

SKU: EB06989

SPECIFICATIONS

Unit Size

100 µg

Synonym /

monoclonal antibody Ki-67|Proliferation-related Ki-67 antigen|Ki-67|KIA|HGNC:7107|antigen identified by

Alias monoclonal antibody Ki-67|MKI67 **Names**

Accession

ID

NP_002408.3

Blocking

EBP06989 **Peptide**

Immunogen

Peptide with sequence CVKKITTRSHRDSED, from the C Terminus of the protein sequence according to

NP 002408.3.

Peptide

Sequence

CVKKITTRSHRDSED

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

Gene ID

4288

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection

Antibody detection limit dilution 1:32000.

Limit

Western

Not yet tested - our routine western blotting protocol does not allow for the detection of proteins >250kDa.

Therefore we cannot recommend an optimal concentration, and the antibody is an aspiring product.

Application

Type

Blot

Pep-ELISA