

GOAT ANTI-KI-67 / MKI67 ANTIBODY

SKU: EB06989



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	monoclonal antibody Ki-67 Proliferation-related Ki-67 antigen Ki-67 KIA HGNC:7107 antigen identified by monoclonal antibody Ki-67 MKI67
Accession ID	NP_002408.3
Blocking Peptide	EBP06989
Immunogen	Peptide with sequence CVKKITTRSHRSED, from the C Terminus of the protein sequence according to NP_002408.3.
Peptide Sequence	CVKKITTRSHRSED
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	4288
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Not yet tested - our routine western blotting protocol does not allow for the detection of proteins >250kDa. Therefore we cannot recommend an optimal concentration, and the antibody is an aspiring product.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

