

GOAT ANTI-KIF26B ANTIBODY

SKU: EB10242



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	KIF26B MGC35030 kinesin family member 26B FLJ10157
Accession ID	NP_060482.2
Blocking Peptide	EBP10242
Immunogen	Peptide with sequence SHYECLSLER, from the internal region of the protein sequence according to NP_060482.2.
Peptide Sequence	SHYECLSLER
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	55083
Mouse Gene ID	269152
Rat Gene ID	305012
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)