

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-KIF26B ANTIBODY**

**SKU:** EB10242

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias KIF26B|MGC35030|kinesin family member 26B|FLJ10157

**Names** 

Accession NP 060482.2 ID

**Blocking** EBP10242

**Peptide** 

Immunogen Peptide with sequence SHYECLSLER, from the internal region of the protein sequence according to NP 060482.2.

**Peptide SHYECLSLER** Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** Refrigerated Instructions

**Predicted** Human, Mouse, Rat

**Species** 

Human

55083 Gene ID

Mouse 269152

**Gene ID Rat Gene ID** 305012

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

**Application** Pep-ELISA **Type**