


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-KINESIN 1 / UKHC, BIOTINYLATED ANTIBODY

SKU: EB05492-B







250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

ubiquitous kinesin heavy chain|kinesin heavy chain|kinesin 1 (110-120kD)|epididymis secretory protein Li **Alias** 61|conventional kinesin heavy chain|UKHC|KNS1|KNS|KINH|HEL-S-61|kinesin family member 5B|KIF5B Names

Accession

NP 004512.1 ID

Blocking

EBP05492-B **Peptide**

Peptide with sequence CQPVAVRGGGGKQV., from the C Terminus of the protein sequence according to **Immunogen**

NP 004512.1.

Peptide Sequence

CQPVAVRGGGGKQV.

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Cow **Species**

Reactive

Human, Dog **Species**

Human

3799 Gene ID

Mouse

16573 Gene ID

Rat Gene ID 117550

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Western

Approx 120kDa band observed in lysates of cell line HeLa (calculated MW of 110kDa according to NP_004512.1). See non-biotinylated parental product's datasheet for further QC data. Recommended

Blot concentration: 1-3µg/ml.

Application

Pep-ELISA, WB **Type**






Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GALLERY IMAGES







250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa