

GOAT ANTI-KLF7 (UKLF) ANTIBODY

SKU: EB05470



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	ubiquitous Kruppel-like transcription factor UKLF Kruppel-like factor 7 (ubiquitous) KLF7
Names	
Accession ID	NP_003700
Immunogen	Peptide with sequence DVLASYISIFQEL-C, from the N Terminus of the protein sequence according to NP_003700.
Peptide Sequence	DVLASYISIFQEL-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Mouse
Human Gene ID	8609
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 35kDa band seen in Mouse Kidney lysate (predicted MW 33kDa). Recommended for use at 2-4µg/ml.