

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-KLK1B27 (MOUSE) ANTIBODY

SKU: EB09189

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias kallikrein 27|glandular kallikrein 27|mGK-27|Klk27|Gk27|kallikrein 1-related peptidase b27|Klk1b27

Names

Accession NP 064664.1 ID

Blocking

EBP09189 **Peptide**

Peptide with sequence C-NDHIPHPEDKSND, from the internal region of the protein sequence according to Immunogen

NP 064664.1.

Peptide

C-NDHIPHPEDKSND Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Mouse

Mouse Gene ID

16619

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:1000.

Limit

Western Preliminary experiments in Mouse Kidney and Testis lysates gave no specific signal but low background (at

Blot antibody concentration up to 1µg/ml).

Application

Pep-ELISA **Type**