

GOAT ANTI-KLK1B27 (MOUSE) ANTIBODY

SKU: EB09189



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	kallikrein 27 glandular kallikrein 27 mGK-27 Klk27 Gk27 kallikrein 1-related peptidase b27 Klk1b27
Accession ID	NP_064664.1
Blocking Peptide	EBP09189
Immunogen	Peptide with sequence C-NDHIPHPEDKSND, from the internal region of the protein sequence according to NP_064664.1.
Peptide Sequence	C-NDHIPHPEDKSND
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Mouse
Mouse Gene ID	16619
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Preliminary experiments in Mouse Kidney and Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)