

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-KLOTHO / KL ANTIBODY

SKU: EB09087

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias KL|klotho

Names

Accession NP 004786.2 ID

Blocking

EBP09087 **Peptide**

Peptide with sequence NQSQVNHTILQYYRC, from the internal region of the protein sequence according to Immunogen

NP 004786.2.

Peptide NQSQVNHTILQYYRC Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

9365 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:1000.

Limit

Preliminary experiments in Human Kidney, Pancreas, Prostate and Testis lysates gave no specific signal but low Western background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field -

Blot have any results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**