


GOAT ANTI-KU70 / G22P1 ANTIBODY

SKU: EB07210



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	thyroid-lupus autoantigen p70 thyroid autoantigen 70kD (Ku antigen) OTTHUMP00000028581 CTC box binding factor 75 kDa subunit ATP-dependent DNA helicase II, 70 kDa
Names	subunit TLAA ML8 KU70 CTCBF CTC75 HGNC:4055 thyroid autoantigen 70kDa (Ku antigen) G22P1
Usage	Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration:
Summary	10ug/ml.
Accession ID	NP_001460.1; NP_001275907.1
Blocking Peptide	EBP07210
Immunogen	Peptide with sequence C-TEKDKN SVNFKN, from the internal region of the protein sequence according to NP_001460.1; NP_001275907.1.
Peptide Sequence	C-TEKDKN SVNFKN
Shipping Instructions	Refrigerated

Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	2547
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, FC

GALLERY IMAGES

