

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-KYNURENINE 3-MONOOXYGENASE (AA29-42) ANTIBODY

SKU: EB12988





Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

KMO|kynurenine 3-monooxygenase (kynurenine 3-hydroxylase)|dJ317G22.1|kynurenine 3-

hydroxylase|kynurenine 3-monooxygenase Names

Accession

ID

NP_003670.2

Blocking

EBP12988 **Peptide**

Peptide with sequence C-KRNFQIDVYEARED, from the internal region (near N terminus) of the protein **Immunogen**

sequence according to NP_003670.2.

Peptide

Sequence

C-KRNFQIDVYEARED

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human

Species Reactive

Human

Species

Human 8564

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:64000.

Limit

Western Approx 55kDa band observed in Human Placenta lysates (calculated MW of 55.8kDa according to

Blot NP_003670.2). Recommended concentration:1-3µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa