



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-KYNURENINE 3-MONOOXYGENASE (AA29-42) ANTIBODY

SKU: EB12988





SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	KMO kynurenine 3-monooxygenase (kynurenine 3-hydroxylase) dJ317G22.1 kynurenine 3-hydroxylase kynurenine 3-monooxygenase
Accession ID	NP_003670.2
Blocking Peptide	EBP12988
Immunogen	Peptide with sequence C-KRNFQIDVYEARED, from the internal region (near N terminus) of the protein sequence according to NP_003670.2.
Peptide Sequence	C-KRNFQIDVYEARED
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	8564
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Approx 55kDa band observed in Human Placenta lysates (calculated MW of 55.8kDa according to NP_003670.2). Recommended concentration:1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GALLERY IMAGES

