

GOAT ANTI-LAG3 ANTIBODY

SKU: EB07683



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	lymphocyte-activation protein 3 CD223 lymphocyte-activation gene 3 LAG-3 LAG3
Names	
Accession ID	NP_002277.3
Blocking Peptide	EBP07683
Immunogen	Peptide with sequence C-QPRVQLDERGRQRGD, from the internal region of the protein sequence according to NP_002277.3.
Peptide Sequence	C-QPRVQLDERGRQRGD
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	3902
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments in Human Lymph Node, Spleen and Tonsil lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA