


Telephone: (650) 697-3600

GOAT ANTI-LAMIN A / LAMIN A/C ANTIBODY

SKU: EB06696



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Charcot-Marie-Tooth disease, axonal, type 2B1|progeria 1 (Hutchinson-Gilford type)|70 kDa lamin|lamin A/C|lamin A\/C|lamin A|LGMD1B|CMT2B1|CMD1A|PRO1|LMNC|LMN1|LDP1|HGPS|FPLD|EMD2|LFP|FPL|LMNA

Names Accession

Alias

ID

NP_733821.1; NP_733822.1

Blocking

Peptide

EBP06696

Immunogen

Peptide with sequence PRTQSPQNCSIM, from the C Terminus of the protein sequence according to NP_733821.1; NP_733822.1.

Product

This antibody is expected to recognise isoforms 1 and 3 (NP_733821.1; NP_733822.1). Comments

Peptide

PRTQSPQNCSIM Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

Predicted

Human, Mouse, Rat

Species Reactive

Human Species

Human **Gene ID** 4000

Product Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$

IHC Results Paraffin embedded Human Prostate. Recommended concentration: 5µg/ml.

ELISA

Antibody detection limit dilution 1:32000. Detection

Limit

Application Pep-ELISA, IHC Type

GALLERY IMAGES

