


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LARP2 / LARP1B ANTIBODY





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias

La ribonucleoprotein domain family member 2|MGC75174|MGC117277|DKFZp686L13217|La ribonucleoprotein

domain family, member 2|LARP2 **Names**

Usage **Summary** Immunofluorescence: Strong expression of the protein seen in the cytoplasm of A431 cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of

HepG2 cells. Recommended concentration: 10ug/ml.

Accession

NP_060548.2

Blocking Peptide

EBP07617









Peptide with sequence C-QGVLHIPKKDLTDE, from the internal region of the protein sequence according to **Immunogen**

NP 060548.2.

Product This antibody is expected to recognise isoform 1 (NP_060548.2) only. **Comments**

Peptide

C-QGVLHIPKKDLTDE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Reactive

Human

Species Human

55132

Gene ID **Product**

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Antibody detection limit dilution 1:128000. Detection

Limit

Application

Pep-ELISA, IF, FC **Type**

GALLERY IMAGES

