


GOAT ANTI-LAVA LAMP (DROSOPHILA MELANOGASTER) ANTIBODY

SKU: EB11446



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names Lva| lvi| lvi| lva-PC| Lva| LVA| Lava-Lamp| lavalamp| lava lamp| lava| Lava| Dmel_CG6450| Dmel\CG6450| CG6450-PC| CG6450| CG#6450| actin-binding protein 2| ABP2

Usage Summary **Immunocytochemistry:** Colabeling with GM130 of Golgi in 4% PFA-fixed Drosophila S2 cells. Recommended concentration: 10ug/ml.

Accession ID NP_525064.1

Blocking Peptide EBP11446

Immunogen	Peptide with sequence C-HTDAQVSAELAKQ, from the internal region of the protein sequence according to NP_525064.1.
Peptide Sequence	C-HTDAQVSAELAKQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Drosophila
Reactive Species	Drosophila
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA, ICC, IF

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

