

GOAT ANTI-LAX1 ANTIBODY

SKU: EB09709



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	linker for activation of x cells OTTHUMP0000034094 LAT-like membrane associated protein LAX FLJ20340 lymphocyte transmembrane adaptor 1 LAX1
Accession ID	NP_060243.2; NP_001129662.1
Blocking Peptide	EBP09709
Immunogen	Peptide with sequence NLTPSAHCINVRAS, from the internal region of the protein sequence according to NP_060243.2; NP_001129662.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_060243.2; NP_001129662.1).
Peptide Sequence	NLTPSAHCINVRAS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	54900
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)