



Telephone: (650) 697-3600



GOAT ANTI-LAX1 ANTIBODY

SKU: EB09709

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

linker for activation of x cells|OTTHUMP0000034094|LAT-like membrane associated **Alias**

protein|LAX|FLJ20340|lymphocyte transmembrane adaptor 1|LAX1 **Names**

Accession

NP_060243.2; NP_001129662.1

Blocking

EBP09709 **Peptide**

Peptide with sequence NLTPSAHCINVRAS, from the internal region of the protein sequence according to Immunogen

NP 060243.2; NP 001129662.1.

Product Comments

This antibody is expected to recognize both reported isoforms (NP_060243.2; NP_001129662.1).

Peptide

NLTPSAHCINVRAS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Species

Human

Human

54900

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:1000.

Limit

Application Pep-ELISA **Type**

> Goat Anti-LAX1 Antibody https://everestbiotech.com/products/goat-anti-lax1-antibody/