



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LCAT ANTIBODY

SKU: EB12992



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa


50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

phospholipid-cholesterol acyltransferase|phosphatidylcholine-sterol acyltransferase|lecithin-cholesterol

acyltransferase|LCAT Names

Accession

ID

NP_000220.1

Blocking

Peptide

EBP12992

Peptide with sequence C-PYDWRLEPGQ, from the internal region of the protein sequence according to **Immunogen**

NP 000220.1.

Peptide

Sequence

C-PYDWRLEPGQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human

Human

Reactive

Species

Human

3931

Gene ID **Product**

Grade

ELISA

Detection

Limit

Antibody detection limit dilution 1:64000.

Western

Approx 55kDa band observed in Human Testis lysates (calculated MW of 49.6kDa according to NP_000220.1). **Blot** This molecular weight is routinely observed by other sources. Recommended concentration: 0.3-1µg/ml.

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa