


GOAT ANTI-LCP1 / PLASTIN-2 (AA277-291) ANTIBODY

SKU: EB12578



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	plastin 2 bA139H14.1 (lymphocyte cytosolic protein 1 (L-plastin)) OTTHUMP00000018362 Lymphocyte cytosolic protein-1 (plasmin) L-plastin (Lymphocyte cytosolic protein 1) (LCP-1) (LC64P) L-plastin lymphocyte cytosolic protein 1 (L-plastin) PLS2 LPL LC64P L-PLASTIN FLJ39956 FLJ26114 FLJ25423 DKFZp781A23186 CP64 LCP1
Accession ID	NP_002289.2
Blocking Peptide	EBP12578
Immunogen	Peptide with sequence HLENAGCNKIGNFST, from the internal region of the protein sequence according to NP_002289.2.
Peptide Sequence	HLENAGCNKIGNFST
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	3936
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 70kDa band observed in lysates of cell lines Daudi and U937 (calculated MW of 70.3kDa according to NP_002289.2). Recommended concentration: 0.5-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

