

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LILRA4 (AA68-78) ANTIBODY

SKU: EB12996



Email: customerservice@vectorlabs.com

Telephone: <u>(650)</u> 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

LILRA4|leukocyte immunoglobulin-like receptor, subfamily A (with TM domain), member 4|CD85g|ILT7|CD85 antigen-like family member G|ILT-7|immunoglobulin-like transcript 7|leukocyte immunoglobulin-like receptor

subfamily A member 4|leukocyte immunog Names

Accession

Alias

NP_036408.3 ID

Blocking

EBP12996 **Peptide**

Immunogen

Peptide with sequence C-SRHILKTLESE, from the internal region of the protein sequence according to

NP 036408.3.

Peptide

C-SRHILKTLESE

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

23547 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Antibody detection limit dilution 1:16000. **Detection**

Limit

Western Approx 55kDa band observed in Human Lymph Node lysates (calculated MW of 55.2kDa according to

Blot NP_036408.3). Recommended concentration: 0.1-0.3µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa