



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-LILRA4 (AA68-78) ANTIBODY

SKU: EB12996





SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	LILRA4 leukocyte immunoglobulin-like receptor, subfamily A (with TM domain), member 4 CD85g ILT7 CD85 antigen-like family member G ILT-7 immunoglobulin-like transcript 7 leukocyte immunoglobulin-like receptor subfamily A member 4 leukocyte immunog
Accession ID	NP_036408.3
Blocking Peptide	EBP12996
Immunogen	Peptide with sequence C-SRHILKTLESE, from the internal region of the protein sequence according to NP_036408.3.
Peptide Sequence	C-SRHILKTLESE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23547
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 55kDa band observed in Human Lymph Node lysates (calculated MW of 55.2kDa according to NP_036408.3). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GALLERY IMAGES

