

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LIPIN 1 (C TERMINUS) ANTIBODY

SKU: EB06663

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias PAP1|KIAA0188|DKFZp781P1796|lipin 1|KIAA0188|LPIN1

Names

Accession NP 663731.1 ID

Blocking

EBP06663 **Peptide**

Immunogen Peptide with sequence C-PFENQDIHSASA, from the C Terminus of the protein sequence according to NP 663731.1.

Peptide C-PFENQDIHSASA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Refrigerated

Instructions

Predicted Human, Mouse, Rat, Dog, Pig, Cow **Species**

Human

23175 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Preliminary experiments in Human Adipose and Lung lysates gave no specific signal but low background (at Western antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any **Blot**

results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**