

GOAT ANTI-LIPIN 1 (C TERMINUS) ANTIBODY

SKU: EB06663



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PAP1 KIAA0188 DKFZp781P1796 lipin 1 KIAA0188 LPIN1
Accession ID	NP_663731.1
Blocking Peptide	EBP06663
Immunogen	Peptide with sequence C-PFENQDIHSASA, from the C Terminus of the protein sequence according to NP_663731.1.
Peptide Sequence	C-PFENQDIHSASA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Human Gene ID	23175
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in Human Adipose and Lung lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA