


Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-LIPIN 2 ANTIBODY**

**SKU:** EB06664



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** lipin 2|KIAA0249|LPIN2

**Names** 

<strong>Immunofluorescence:</strong> Strong expression of the protein seen in the cytoplasm of A431 and Usage U2OS cells. Recommended concentration: 10µg/ml. <strong>Flow Cytometry:</strong> Flow cytometric **Summary** 

analysis of A431 cells. Recommended concentration: 10ug/ml.

**Accession** 

NP 055461.1

**Blocking Peptide** 

ID

EBP06664









Peptide with sequence CYWRDPIPEVDLDDLS, from the C Terminus of the protein sequence according to **Immunogen** NP 055461.1.

**Peptide CYWRDPIPEVDLDDLS Sequence** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Pig **Species** 

Reactive

Human **Species** 

Human

9663

Gene ID **Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:8000.

Limit

Application

Pep-ELISA, IF, FC **Type** 

## **GALLERY IMAGES**

