

GOAT ANTI-LIPOCALIN 2 / NGAL (MOUSE) ANTIBODY

SKU: EB08266

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	NGAL AW212229 24p3 RP23-161B9.11 NGAL lipocalin 2 Lcn2
Accession ID	NP_032517.1
Blocking Peptide	EBP08266
Immunogen	Peptide with sequence C-HRYPQVQSYNVQ , from the internal region of the protein sequence according to NP_032517.1.
Peptide Sequence	C-HRYPQVQSYNVQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Mouse
Reactive Species	Mouse
Mouse Gene ID	16819
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 23kDa band observed in Mouse Kidney lysates (calculated MW of 22.6kDa according to NP_005555.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

