

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-LIPOCALIN 2 / NGAL (MOUSE) ANTIBODY**

**SKU:** EB08266



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa 50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

NGAL|AW212229|24p3|RP23-161B9.11|NGAL|lipocalin 2|Lcn2

Names

Accession

NP\_032517.1 ID

**Blocking** 

EBP08266 **Peptide** 

Peptide with sequence C-HRYPQVQSYNVQ, from the internal region of the protein sequence according to **Immunogen** 

NP 032517.1.

**Peptide** Sequence

C-HRYPQVQSYNVQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Mouse **Species** 

Reactive

Mouse **Species** 

Mouse

16819 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Grade **ELISA** 

Antibody detection limit dilution 1:8000. **Detection** 

Limit

Western Approx 23kDa band observed in Mouse Kidney lysates (calculated MW of 22.6kDa according to NP\_005555.2).

**Blot** Recommended concentration: 1-3µg/ml.

**Application** 

Pep-ELISA, WB Type

## **GALLERY IMAGES**



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa