

# GOAT ANTI-LIPOPHILIN B / SCGB1D2 ANTIBODY

SKU: EB06809



## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	SCGB1D2 secretoglobin, family 1D, member 2 HGNC:18396 LIPB LPHB lipophilin B lipophilin B (uteroglobin family member), prostatein-like prostatein-like lipophilin B
<b>Accession ID</b>	NP_006542.1
<b>Blocking Peptide</b>	EBP06809
<b>Immunogen</b>	Peptide with sequence CTDQMSLQKRSLIAE, from the internal region of the protein sequence according to NP_006542.1.
<b>Peptide Sequence</b>	CTDQMSLQKRSLIAE
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Human Gene ID</b>	10647
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:32000.
<b>Western Blot</b>	Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 9.9kDa according to NP_006542.1. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)