

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LIPOPHILIN B / SCGB1D2 ANTIBODY

SKU: EB06809

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / SCGB1D2|secretoglobin, family 1D, member 2|HGNC:18396|LIPB|LPHB|lipophilin B|lipophilin B (uteroglobin family

Alias member), prostatein-like|prostatein-like lipophilin B

Names

Accession NP 006542.1 ID

Blocking

EBP06809 Peptide

Peptide with sequence CTDQMSLQKRSLIAE, from the internal region of the protein sequence according to Immunogen

NP 006542.1.

Peptide

CTDQMSLQKRSLIAE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

10647 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:32000.

Limit

Blot

Western

Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 9.9kDa according to NP_006542.1. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been

reported?

Application

Pep-ELISA

Type





Telephone: (650) 697-3600

