



GOAT ANTI-LIPOPHILIN B / SCGB1D2 ANTIBODY

SKU: EB06809



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SCGB1D2 secretoglobin, family 1D, member 2 HGNC:18396 LIPB LPHB lipophilin B lipophilin B (uteroglobin family member), prostatein-like prostatein-like lipophilin B
Accession ID	NP_006542.1
Blocking Peptide	EBP06809
Immunogen	Peptide with sequence CTDQMSLQKRSLIAE, from the internal region of the protein sequence according to NP_006542.1.
Peptide Sequence	CTDQMSLQKRSLIAE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	10647
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 9.9kDa according to NP_006542.1. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
Application Type	Pep-ELISA



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Data Sheet](#)