


GOAT ANTI-LIS1 / PAFAH1B1 ANTIBODY

SKU: EB05439



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MDS LIS2 platelet-activating factor acetylhydrolase, isoform Ib, subunit 1 (45kDa) platelet-activating factor acetylhydrolase, isoform Ib, alpha subunit (45kD) Platelet-activating factor acetylhydrolase, isoform 1B, alpha subunit lissencephaly 1 protein PAFAH MDCR platelet-activating factor acetylhydrolase, isoform Ib, alpha subunit 45kDa LIS1 PAFAH1B1
Accession ID	NP_000421.1
Blocking Peptide	EBP05439
Immunogen	Peptide with sequence TGSVDQTVKWECD, from the C Terminus of the protein sequence according to NP_000421.1.
Peptide Sequence	TGSVDQTVKWECD
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Rat
Reactive Species	Human, Rat
Human Gene ID	5048
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Brain (Cortex), Thyroid and Pancreas. Recommended concentration: 3.75µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 40kDa band observed in Human Lung and Brain (Cerebellum) and in Rat Ovary lysates (calculated MW of 46.6kDa according to Human NP_000421.1 and to Rat NP_113951.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

