



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LMO6 ANTIBODY

SKU: EB08212



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias prickle homolog 3 (Drosophila)|PRICKLE3|LMO6|LIM domain only 6

Names

Accession

NP_006141.2 ID

Blocking

EBP08212 **Peptide**

Peptide with sequence HPSRRRHYQCD, from the internal region of the protein sequence according to **Immunogen**

NP 006141.2.

Peptide

HPSRRRHYQCD Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

4007 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western Approx 70kDa band observed in lysates of cell line Daudi (calculated MW of 68.6kDa according to

Blot NP_006141.2). Recommended concentration: 0.1-0.3µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: <u>(650)</u> 697-3600

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa