




GOAT ANTI-LMX1B ANTIBODY

SKU: EB08844



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	LIM homeo box transcription factor 1, beta NPS1 MGC142051 MGC138325 LMX1.2 LIM homeobox transcription factor 1, beta LMX1B
Accession ID	NP_002307.1
Blocking Peptide	EBP08844
Immunogen	Peptide with sequence QQALTTCYFRDRK , from the internal region of the protein sequence according to NP_002307.1.
Peptide Sequence	QQALTTCYFRDRK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated



Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	4010
Mouse Gene ID	16917
Rat Gene ID	114501
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Testis. Recommended concentration: 2.5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

