

GOAT ANTI-LNP1 ANTIBODY

SKU: EB09643



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	leukemia NUP98 fusion partner 1 LNP1
Accession ID	NP_001078920.1
Blocking Peptide	EBP09643
Immunogen	Peptide with sequence SSDCHPRRSHEDQ, from the internal region of the protein sequence according to NP_001078920.1.
Peptide Sequence	SSDCHPRRSHEDQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Cow
Human Gene ID	348801
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)