Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-LNP1 ANTIBODY

SKU: EB09643

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias leukemia NUP98 fusion partner 1|LNP1

Names

Accession NP 001078920.1 ID

Blocking

EBP09643 **Peptide**

Peptide with sequence SSDCHPRRHSHEDQ, from the internal region of the protein sequence according to Immunogen

NP 001078920.1.

Peptide SSDCHPRRHSHEDQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Dog, Cow **Species**

Human

348801 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:128000.

Limit

Application Pep-ELISA **Type**