

GOAT ANTI-LNX ANTIBODY

SKU: EB05962



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	multi-PDZ-domain-containing protein MPDZ ligand of numb-protein X LNX
Accession ID	NP_116011
Blocking Peptide	EBP05962
Immunogen	Peptide with sequence C-RITLTIVSWPGTFL, from the C Terminus of the protein sequence according to NP_116011.
Peptide Sequence	C-RITLTIVSWPGTFL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	84708
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Human Kidney extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA